


Uni-Lock System





Project Overview

Scope


- Provide a unique injector system with Prima Tech Injectors for use with injectable products that offers reliability, consistency and ease of use to customers

Uni-Lock Draw off

- Accommodates 20mm, 30mm and 33mm Bottles
- Eliminates typical draw off leakage due to loose fit on bottles
- Allows for smooth change out of bottles versus all other systems available today
- Fast Fit offers Built in Bottle Protection when hanging on chute side




Uni-Lock Draw off with Collar





Uni-Lock Draw off with Fast Fit





Bottle Set Ups



100ml



200ml